en

SURGICEL™ ORIGINAL and SURGICEL™ NU-KNIT

ABSORBABLE HAEMOSTATS

(OXIDISED REGENERATED CELLULOSE)

Description

SURGICEL™ Absorbable Haemostat is a sterile absorbable knitted fabric prepared by the controlled oxidation of regenerated cellulose. SURGICEL™ Absorbable Haemostat is 100% composed of oxidised regenerated cellulose. The fabric is white with a pale-yellow cast and has a faint, caramel-like aroma. It is strong and can be sutured or cut without fraying. A slight discoloration may occur with age, but this does not affect performance. SURGICEL™ Absorbable Haemostat is intended for use only by healthcare professionals who are trained in the surgical procedures and techniques requiring the use of an Absorbable Haemostat.

The clinical benefit to be expected is the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective.

A summary of safety and clinical performance can be found at the following link (upon activation): https://ec.europa.eu/tools/eudamed.

Indications / Intended Use

SURGICEL™ Haemostat is used adjunctively in surgical procedures to assist in the control of capillary, venous, and small arterial haemorrhage when ligation or other conventional methods of control are impractical or ineffective.

SURGICEL™ Haemostat can be cut to size for use in endoscopic procedures (see Figures 1, 2A, 2B and 3).

Figure 1. SURGICEL™ *Haemostat should be cut to the appropriate size for endoscopic placement. Standard endoscopic procedures should be used up to the point of placement of the absorbable haemostat. Grasp the* SURGICEL™ *Haemostat at one corner.*

Figures 2A and 2B. *Slowly push the grasping instrument and material into the cavity.*

Figure 3. *With the use of grasping instruments in a second and/or third auxiliary site, placement can be made, and the material repositioned as needed.*

Contraindications

- Although packing or wadding sometimes is medically necessary, SURGICEL™ Haemostat should not be used in this manner, unless it is to be removed after haemostasis is achieved.
- SURGICEL™ Haemostat should not be used for implantation in bone defects, such as fractures, since there is a possibility of interference with callus formation and a theoretical chance of cyst formation.

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