

FICHA TÉCNICA

PANTALLA DE PROTECCIÓN FACIAL BESTDENT		
REFERENCIAS	78433: Máscara de protección + montura (Modelo R66). 78434: Reposición máscara protección.	
FABRICANTE/ RESPONSABLE	LARIDENT SRL Via Lamaneigra 12R - 16030 Tribogna (GE) ITALIA +39 0185 939481	
CLASIFICACIÓN	Equipo de Protección individual (EPI) _ Clase II, de acuerdo la EN 166:2024 y el Reglamento 2016/425.	
DESCRIPCIÓN	- Pantalla protectora facial de policarbonato con marco de clip. - El visor cubre nariz y boca para una protección completa. - Posee una doble capa protectora y es muy transparente, garantizando una excelente visibilidad.	
COMPOSICIÓN	Policarbonato.	
ENVASE	78433: 1 montura + 12 máscaras 78434: 12 máscaras.	
FINALIDAD PREVISTA	Protección facial.	
INSTRUCCIONES DE USO	Se suministra en condición "NO ESTÉRIL". Autoclavable. Marco esterilizable en autoclave hasta 121°C. Instrucciones de montaje disponibles dentro de la caja de cartón. El visor protector NO es resistente a soluciones alcohólicas. ADVERTENCIA: Solo para uso como pantalla facial contra salpicaduras de líquidos acuosos (agua y saliva).	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

FACE SHIELD PROTECTION BESTDENT		
REFERENCES	78433: Face Shield protection + monture (Model R66). 78434: Face Shield replacements.	
MANUFACTURER/ RESPONSIBLE	LARIDENT SRL Via Lamaneigra 12R - 16030 Tribogna (GE) ITALY +39 0185 939481	
CLASSIFICATION	Personal Protective Equipment (PPE) _ Class II, according to EN 166:2024 and Regulation 2016/425.	
DESCRIPTION	- Polycarbonate face shield with clip-on frame. - The visor covers nose and mouth for complete protection. - It has a double protective layer and is highly transparent, ensuring excellent visibility.	
COMPOSITION	Polycarbonate.	
PACKAGE	78433: 1 monture + 12 masks. 78434: 12 masks.	
INTENDED USE	Face protection.	
ISTRUCCIONS FOR USE	Supplied in “NON-STERILE” condition. Autoclavable. Frame autoclavable up to 121°C. Mounting instructions available inside the carton box. The face shield is NOT resistant to alcoholic solutions. WARNING: Only for use as a face shield against splashes of aqueous liquids (water and saliva).	
STORAGE AND TRANSPORT	It does not require special storage and transportation conditions.	