

FICHA TÉCNICA

CEMENTO IONÓMERO DE VIDRIO PROCLINIC EXPERT		
REFERENCIAS	78555	
FABRICANTE / RESPONSABLE	P.L. Superior Dental Materials GmbH Stiller Weg 15a 22607-Hamburg Web: www.sc-polymer.com	
CLASIFICACIÓN	Producto Sanitario Clase IIa, reglamento (UE) 2017/745 MDR.	
DESCRIPCIÓN	- Cemento de ionómero de vidrio clásico radiopaco pensado para cementar una amplia variedad de materiales protésicos. - Ofrece una buena adhesión química y una gran biocompatibilidad.	
COMPOSICIÓN	- Polvo: vidrio de silicato de estroncio fluoroaluminio, ácidos poliacrílicos secos, pigmentos. - Líquido: solución acuosa de ácidos poliacrílicos.	
PRESENTACIÓN	15 g polvo + 10 ml líquido.	
VIDA ÚTIL	2 años.	
FINALIDAD PREVISTA	- Cementación de coronas, puentes, inlays y onlays, postes endodónticos o bandas de ortodoncia fabricados en todo tipo de metal y cerámica. - Como base cavitaria en restauraciones de composite y amalgama.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	Conserve Proclinic Expert Glass Ionomer Cement en un lugar fresco y oscuro entre 4 y 25°C (39-78 °F). La temperatura no debe superar los 25 °C (78 °F). Proteja de la humedad.	

TECHNICAL DATA SHEET

GLASS IONOMER CEMENT PROCLINIC EXPERT		
REFERENCES	78555	
MANUFACTURER	P.L. Superior Dental Materials GmbH Stiller Weg 15a 22607-Hamburg Web: www.sc-polymer.com	
CLASSIFICATION	Class IIa Medical Device, regulation (EU) 2017/745 MDR.	
DESCRIPTION	- Classic radiopaque glass ionomer cement designed to cement a wide variety of prosthetic materials. - It offers good chemical adhesion and great biocompatibility.	
COMPOSITION	- Powder: strontium fluoro alumino silicate glass, dried polyacrylic acids, pigments. - Liquid: aqueous solution of polyacrylic acids	
PACKAGE	15 g powder + 10 ml liquid.	
SHELF LIFE	2 years.	
INTENDED USE	- Cementation of crowns, bridges, inlays and onlays, endodontic posts or orthodontic bands made of all types of metal and ceramic. - As a cavity base in composite and amalgam restorations.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	Store Proclinic Expert Glass Ionomer Cement in a cool and dark place at 4-25 °C (39-78 °F). Temperature should not exceed 25 °C (78 °F).	