

FICHA TÉCNICA

BRUÑIDOR 27/29 BESTDENT		
REFERENCIA	80245	
FABRICANTE / RESPONSABLE	MARKONE SURGICAL Ugoki Road, Adalat Garah, Sialkot 51330, Pakistan	
CLASIFICACIÓN	Producto Sanitario Clase I, Reglamento MDR (EU) 2017/745.	
DESCRIPCIÓN	Bruñidor 27/29 <u>Características:</u> Instrumento con doble parte activa terminación en bola y en pera. No presentan alteraciones con los desinfectantes químicos. Dureza mínima entre 40-44hrc	
ESPECIFICACIÓN	<u>Tamaños de las bolas:</u> 1,6 mm / 3 mm	
COMPOSICIÓN	Acero inoxidable quirúrgico.	
PRESENTACIÓN	Contiene una unidad.	
VIDA ÚTIL	La vida útil del instrumento depende de la frecuencia de uso, el cuidado del usuario y los métodos de reprocesamiento adecuados.	
ESTERILIZACIÓN	Todos los instrumentos se deben examinar, limpiar y esterilizar antes de cada uso. Esterilizar a 177 °C en autoclave durante al menos 4 minutos, seguido de un secado adecuado.	
FINALIDAD PREVISTA	Instrumento indicado para la condensación de las obturaciones hechas en composite.	
ALMACENAMIENTO Y TRANSPORTE	No requiere condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

27/29 BESTDENT POLISHER		
REFERENCE	80245	
MANUFACTURER	MARKONE SURGICAL Ugoki Road, Adalat Garah, Sialkot 51330, Pakistan	
CLASSIFICATION	Class I Medical Device, MDR Regulation (EU) 2017/745.	
DESCRIPTION	Burnisher 27/29 <u>Features:</u> Instrument with double active part, ball and pear-shaped tip. Not affected by chemical disinfectants. Minimum hardness between 40-44 HRC.	
SPECIFICATION	<u>Ball sizes:</u> <u>1.6 mm / 3 mm</u>	
COMPOSITION	Surgical stainless steel.	
PACKAGE	Contains one unit.	
SHELF LIFE	The service life of the instrument depends on the frequency of use, user care, and proper reprocessing methods.	
STERILIZATION	All instruments should be examined, cleaned, and sterilized before each use. Sterilize at 177°C in an autoclave for at least 4 minutes, followed by adequate drying.	
INTENDED USE	Instrument indicated for the condensation of composite fillings.	
STORAGE AND TRANSPORT	No special storage or transportation conditions are required	