

FICHA TÉCNICA

POSTES FIBRA DE CARBONO REPOSICIÓN		
REFERENCIAS	7796: 1,0 MM. 7797: 1,2 MM. 7798: 1,4 MM. 7799: 1,6 MM.	
FABRICANTE / RESPONSABLE	BIOLOREN S.R.L Via Alessandro Volta, 59 21047 Saronno Italy info@bioloren.com	
CLASIFICACIÓN DEL PRODUCTO	Producto Sanitario Clase IIa, Directiva 93/42/CEE.	
DESCRIPCIÓN	Fibra de Carbono que representan una importante innovación para la reconstrucción protésica del diente desvitalizado. Características: Módulo de elasticidad similar a la dentina, ausencia de fractura radicular, estructura con fibra longitudinal, ausencia de oxidación y corrosión, forma cilindro-cónica. Ventajas: Fácil remoción del poste en caso de retratamiento, estabilidad y biocompatibilidad, ausencia de fracturas en el poste, uniforme distribución del estrés oclusal.	
ESPECIFICACIÓN	<u>Medidas:</u> 7796: 1,0 MM. 7797: 1,2 MM. 7798: 1,4 MM. 7799: 1,6 MM.	
COMPOSICIÓN	Fibra de carbono y matriz resinosa epóxida.	
PRESENTACIÓN	8 unidades por envase.	
FINALIDAD	Están diseñados para ser utilizados en tratamientos odontológicos, específicamente en la reconstrucción de dientes endodonciados.	
INSTRUCCIONES DE USO	Ver instrucciones en el producto.	
ALMACENAMIENTO Y TRANSPORTE	No se requieren condiciones especiales de almacenamiento y transporte.	

TECHNICAL DATA SHEET

PROCLINIC CARBON FIBER POSTS REPLACEMENT		
REFERENCES	7796: 1,0 MM. 7797: 1,2 MM. 7798: 1,4 MM. 7799: 1,6 MM.	
MANUFACTURER	BIOLOREN S.R.L Via Alessandro Volta, 59 21047 Saronno Italy info@bioloren.com	
CLASSIFICATION OF THE PRODUCT	Class IIa Medical Device, according to Directive 93/42/EEC.	
DESCRIPTION	Carbon Fiber that represent an important innovation for the prosthetic reconstruction of devitalized teeth. Features: Modulus of elasticity similar to dentin, absence of root fracture, structure with longitudinal fiber, absence of oxidation and corrosion, cylindrical-conical shape. Advantages: Easy removal of the post in case of retreatment, stability and biocompatibility, absence of fractures in the post, uniform distribution of occlusal stress.	
SPECIFICATION	<u>Measures:</u> 7796: 1,0 MM. 7797: 1,2 MM. 7798: 1,4 MM. 7799: 1,6 MM.	
COMPOSITION	Carbon fiber and epoxy resin matrix.	
PACKAGING	8 units per container.	
INTENDED USE	They are designed to be used in dental treatments, specifically in the reconstruction of endodontic teeth.	
INSTRUCTIONS FOR USE	See instructions on the product.	
HANDLING AND STORAGE CONDITIONS	No special storage and transport conditions are required.	